ANEX I: CE CERTIFICATE

AUXILAB S.L.

CE DECLARATION OF CONFORMITY

HORIZONTAL ELECTROFORESIS CELL of Auxilab,S.L for the Directive of machines (89/392/CEE modified) and the regulations adopted for its transposition

NAME OF THE MANUFACTURER / IMPORTER:

AUXILAB, S.L.

ADDRESS:

Polígono Morea Norte, 8
31191 Beriáin (Navarra)

WE STATE THAT:

HORIZONTAL ELECTROFORESIS CELL
53001001

Are designed and manufactured according to:

- ◆ Directive 89/392/CEE, including the modifications and the national regulations that transpose them.
- ◆ Directive 73/23/CEE, modified over the electric security.
- ◆ Directive 89/336/CEE, modified over the electromagnetic compatibility.
- ◆ And that the following harmonized rules have applied (or part of them):
UNE 292-1/-2/-2/A1, UNE-EN 1050, UNE-EN 614-1, UNE-EN 1037, UNE-EN 1088, UNE-EN 547, UNE-EN 953, UNE-EN 294, UNE-EN 418, UNE-EN 894-1, UNE-EN 894-2, UNE-EN 954-1, UNE-EN 60204-1, UNE 61010-1/A2, UNE-EN61010-2-051.

BERIAIN 04th JUNE 2007


Signed by: ALFONSO AINCIBURU SANZ
DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071
Internet: www.auxilab.es · Email: correo@auxilab.es

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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INDICE DE IDIOMAS

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7. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the vertical electrophoresis cell it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

"

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ Never autoclave or dry-heat sterilize the apparatus or components.
- ◆ Do not expose the apparatus or components to phenol, benzene, acetone, halogenated hydrocarbon solvents or undiluted laboratory alcohols.
- ◆ Avoid prolonged exposure of the apparatus or components to the UV light.
- ◆ Periodically make a simple and routine inspection of the electrical components to ensure both the safety and good performance of the equipment:
 - If the banana connectors of either the gel support or the power cords can be rocked easily or present signs of corrosion, replace them.
 - If power cords show any sign of wear and damage replace them immediately.

Cleaning

- ◆ After each use, wash all components gently with water and non abrasive detergents, and rinse well in deionized water. **Note: be careful not to damage the electrode wire when cleaning the gel support.**
- ◆ Wipe dry with a soft cloth or allow to air dry.
- ◆ Never use scourers or substances that can grate as they damage the electrophoresis cell and produce an early ageing of the equipment.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

6. TROUBLESHOOTING

PROBLEM	COMMENTS
The glass plates crack	<ul style="list-style-type: none"> - The wedge to hold the gel is too tight. - The voltage or current setting may be too high, causing overheating. - Electrophoresis buffer may be too concentrated, leading to a too high current and overheating at a constant voltage.
The gel dye front is not straight	<ul style="list-style-type: none"> - Verify that the cell is level during electrophoresis - The gel may be overheated. Reduce the voltage or current. Use chilled electrophoresis buffer.
Bands are distorted	<ul style="list-style-type: none"> - Samples may contain an excess of salt. Dialyze or desalt before loading on gel. - Samples may be too concentrated. Dilute samples or reduce voltage. - Samples may contain precipitate material. Centrifuge or filter before loading on gel. - The upper surface of the resolving gel may not be flat. Use water-saturated n-butanol to overlay the gel solution, rather than water or buffer. - Sample wells may contain small fragments of polyacrylamide. Rinse wells thoroughly before loading samples.
Lanes are broader at the bottom than the top	<ul style="list-style-type: none"> - This is normal if adjoining lanes are not loaded with similar samples. It is also normal for polyacrylamide gradient gels.
Dye migration is slower than expected	<ul style="list-style-type: none"> - The electrophoresis buffer may be too concentrated. Check buffer preparation procedure. If the concentration is high, at constant current the voltage will be lower than usual. - Check voltage or current settings on power supply. - Check for secure connections of lid and power cords.
Dye migration is slower than expected	<ul style="list-style-type: none"> - The electrophoresis buffer may be too diluted. Check buffer preparation procedure. - Check voltage or current settings on power supply.

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


NOTE: acrylamide is a powerful neurotoxin. Always wear gloves when handling solutions, gels or equipment that has been in contact with this substance and strictly observe the safety regulations indicated for its manipulation and disposal of wastes.

1. Resolving gel

CHART 3. To prepare 10 mL of acrylamide solutions of different percentages:

Component	Volume (mL) according to acrylamide %			
	7.5 %	10 %	12.5 %	15 %
Acrylamide stock solution [29% (w/v) acrylamide, 1% (w/v) bisacrylamide in deionized water]	2.5	3.3	4.2	4.0
1.5 M Tris-HCl (pH 8.8)	2.5	2.5	2.5	2.5
10% (w/v) SDS	0.1	0.1	0.1	0.1
Deionized water	4.8	4.0	3.1	2.3
10% (w/v) ammonium persulphate	0.1	0.1	0.1	0.1
TEMED	0.01	0.01	0.01	0.01

2. Stacking gel

CHART 4. To prepare 10 mL of 4 % acrylamide solution:

Component	Volume (mL)
Acrylamide stock solution [29% (w/v) acrylamide, 1% (w/v) bisacrylamide in deionized water]	1.3
1.5 M Tris-HCl (pH 8.8)	2.5
10% (w/v) SDS	0.1
Deionized water	6.0
10% (w/v) ammonium persulphate	0.1
TEMED	0.001

3. Electrophoresis buffer

CHART 5. Final concentration of the different components of the electrophoresis buffer. Final pH must be 8.3.

Component	Final concentration
Tris base	24.8 M
Glycine	192 mM
SDS	0.1% (w/v)

red to one of the glass plates.

- Transfer the gel to a proper container for further staining or electroblotting.

Note: for electroblotting, please refer to Nahita Transfer Module (ref. 53001501) user's manual.

- Clean the tank and all the accessories with deionized water. Allow to air dry.

Security

- ♦ The vertical electrophoresis cell must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ♦ Put the electrophoresis cell on top of a horizontal, plane and stable table.
- ♦ Do not put the electrophoresis cell near any warm supply (burners, blowlamps...). Avoid vibrations.
- ♦ Vertical electrophoresis cell is intended to use with a DC power supply. An inadequate manipulation may cause burns, hurts and damages to the electrical installation.
- ♦ Never connect the electrophoresis cell to the power supply before starting the experiment. The electrical current supplied by the power supply is transmitted to the cell through the leads and the lid assembly.
- ♦ Never use the electrophoresis cell without the lid.
- ♦ Always turn off the power supply before opening and removing the lid.
- ♦ Never operate damaged or leaking equipment.
- ♦ Certain reagents indicated for use in this manual are of a hazardous nature. User should strictly observe the safety regulations indicated for such products.
- ♦ The equipments used in these procedures (e.g. power supplies) should be used following the manufacturer safety recommendations.

5. FORMULATION FOR GELS AND ELECTROPHORESIS BUFFER

Polyacrylamide gels are generated by polymerization of acrylamide as a effect of a cross-linking agent, the bis-acrylamide, and in the presence of an initiator and a catalytic agent. Usually, TEMED (N,N,N,N'-tetramethylethylenediamine) is used as the initiator and the persulphate ion ($S_2O_8^{2-}$) added in the form of ammonium persulphate as the catalytic agent.

Gel porosity depends on the relative proportions of polyacrylamide and bis-acrylamide; a higher bis-acrylamide vs. acrylamide proportion generates a smaller pore. Usually, gels are designated with the acrylamide/bisacrylamide % that contain, thus, this value determinates the separation range of the gel. Most of proteins are well separated in 5-10 % gels; a lower percentage (larger pore) is better to separate larger proteins.

Protein electrophoresis is usually performed in denaturalizing conditions, it means, proteins lose their three-dimensional structure and migrate completely denaturalized. Migration is therefore, proportional to the charge and size of the molecule but not to its shape. The most commonly used denaturalizing agent is the sodium dodecyl sulphate (SDS).

The formulas for preparation of electrophoresis buffer and discontinuous denaturing gels for protein separation are shown below.

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
must have enough glycerol to be denser than the electrophoresis buffer.

Well visualization may be difficult; to facilitate localization of wells you can use two methods:

- a. Before putting the gel support into the tank and with the comb still placed between both glass plates, outline the wells with a water resistant marking pen on the non-notched glass plate.
- b. Illuminate the bench area with a light source placed directly behind the electrophoresis cell.

CHART 2 Maximum volume of sample depending on the type of comb and the gel thickness.

<i>Gel thickness (mm)</i>	<i>Nº teeth</i>	<i>Vol/well (µL)</i>
1.0	10	50
	15	30
1.5	10	75
	15	45

Note: volumes given are approximate.

Note: To obtain the best results and promote a uniform running of the stacking front, load sample buffer in all the wells whether or not they contain samples

- Fill the tank with approximately 250 mL of electrophoresis buffer. Buffer level must not exceed the wells.

♦ Electrophoresis

- Cover the electrophoresis tank with the lid. ***Note: The lid is especially designed to fit in a unique position thus assuring that the black lead is always connected to the black electrode and the red lead to the red electrode.***

- Plug the power leads to a power supply.

- Switch the power supply on and set the desired parameters. Standard SDS gel electrophoresis is conveniently performed 125 to 200 V (constant voltage) for 45 min to 1 h. ***Note: The first time the apparatus is used for particular type of gel, monitor the progress of electrophoresis frequently.***

- When electrophoresis is complete (the marker dye front has reached the bottom of the gel), turn the power supply off and disconnect the leads.

♦ After electrophoresis

- Remove the lid

- Gently pull out the wedge that holds the gel; it may be slippery in detergent-containing solutions. Take the gel support out from the electrophoresis tank.

- Properly discard the electrophoresis buffer. Do not reuse the buffer.

- Take the glass plate assembly with the gel out from the gel support and put them on a paper towel.

- Use a thin spatula or similar to separate the glass plates. The gel will be adhe

Gel thickness	Acrylamide volume			
	1 gel		2 gels	
	Discontinuous	Continuous	Discontinuous	Continuous
1.0 mm	3.4 mL	6.7 mL	6.7 mL	13.5 mL
1.5 mm	5 mL	10 mL	10 mL	20 mL

- Carefully pipet the resolving acrylamide gel filling the space between both glass plates until 1 or 2 mm above the desired level. Overlay the solution with a few drops of water-saturated n-butanol to eliminate bubbles and separate the gel from the air oxygen (inhibitor of polymerization).

- Wait until the gel polymerizes (15-30 min); polymerization is faster at 4 °C. Note: keep the assembly vertical while the acrylamide polymerizes.

- After the resolving gel has polymerized, pour off any unpolymerized solution and n-butanol. Rinse the top of the gel with deionized water and drain thoroughly.

- Add the stacking gel solution to within 2 mm of the top of the notched plate. Eliminate any bubble formed.

- Insert completely the comb. Be careful not to trap bubbles under and between the comb teeth.

- Allow the acrylamide solution to polymerize completely.

2. Continuous buffer system gels

- Prepare the volume of acrylamide solution depending on gel thickness (Chart 1).

- Carefully pipet the acrylamide solution filling the space between both glass plates to within 2 mm of the top of the notched plate. Eliminate any bubble formed.

- Insert completely the comb. Be careful not to trap bubbles under and between the comb teeth.

- Allow the acrylamide solution to polymerize completely.

◆ Assembly for gel electrophoresis and sample loading

Gel electrophoresis is performed on the same support in which it has been cast, minimizing gel manipulation and speeding up and making work easier.

Proceed as follows:

- Turn the handles of the gel casting device to their original position.

- Take off the support with the polymerized gels from the gel casting device and put it into the electrophoresis tank. **Note: the special design of the tank and the gel support allows to fit them at the correct position.**

- Carefully remove the comb by pulling it vertically and rinse with deionized water.

- Fill the central part of the gel support with approximately 150 mL of electrophoresis buffer. Make sure that the wells are completely covered with buffer.

- Carefully rinse the wells with deionized water or electrophoresis buffer to eliminate bubbles or acrylamide fragments. Use a syringe or pipette.

- Load the samples into the wells with a micropipette. The sample volume that can be loaded in each well depending on gel thickness is shown in chart 2. **Note: samples**

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
- Take the exclusive gel casting device and put on its base a silicon rubber cushion. **Note: The cushion must be perfectly clean and dry. After gel casting, remove the silicon rubber cushion, rinse it with water and allow to air dry. Do not use any other method to clean the silicon cushion.**

- Put the gel support into the gel casting device. The base of the gel support must be perfectly sealed with the silicon rubber cushion to avoid acrylamide solution leaks

- Press inwards the white handles placed on each side of the gel casting device; turn them simultaneously until the arrows point at the correct position depending on the gel thickness (each handle presents marks with 1.5 or 1.0 indications that refer to the gel thickness; if a 1.5 mm thick gel is going to be cast, handles should be turned until the arrows point at the 1.5 marks).

◆ Gel casting

Note: Do not move the handles of the gel casting device during gel casting process since the assembly could be loosened resulting in acrylamide leaks.

1. Discontinuous buffer system gels

This system is specially recommended to analyze diluted samples without resolution lost. Gel casting must be performed with two different buffers. The first buffer (stacking gel) assures the accumulation of all proteins loaded in each well at the migration front. When the migration front reaches the second buffer (resolving gel), the separation of the different proteins begins. The first gel, stacking gel, presents a larger pore (lower acrylamide/bisacrylamide percentage) and presents a more acid pH than the resolving gel.

For casting this type of gels proceed as follows:

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CHART 1. Volume of acrylamide solution that must be prepared depending on gel thickness. In discontinuous gels, the indicated volume refers to the volume of solution that must be prepared for each type of gel (stacking or resolving gel).

- 1 thick glass plate for single gel electrophoresis
- Set of combs:
 - 2 combs 1.0 mm thick and with 10 teeth
 - 2 combs 1.0 mm thick and with 15 teeth
 - 2 combs 1.5 mm thick and with 10 teeth
 - 2 combs 1.5 mm thick and with 15 teeth
- Warranty
- User's manual

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

4.2 Operating mode


Note: the equipment is exclusively designed for laboratory purposes; ambient temperature must be 0-40 °C and relative humidity must not exceed the 80%. Prior to use, put the electrophoresis cell on a plane, smooth and vibration-free surface.

◆ Assembly for gel casting

- Make sure that the flat glass plates, the notched glass plates, the combs, the wedges and the gel support are perfectly clean and dry before assembling the different components.

- Put a notched glass plate on a flat plate, so as the fixed spacers stay between both plates. The edges of both plates must be perfectly aligned; to get this, put the plates vertically on a plane surface and make sure that the lateral edges are aligned.

- Put both plates into a compartment of the gel support so as the notched plate faces inside. Holding the plates, insert a plastic wedge (the flan side is in contact with the glass plates) and push it downwards to perfectly hold the glass plates.

Note: if a single gel electrophoresis is going to be performed, put the thick glass plate into the other compartment of the gel support and hold it with the other plastic wedge. In case of a double gel electrophoresis, proceed as previously indicated to prepare the second gel.

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
- **Combs:** the equipment is supplied with a set of 8 combs of 1.0 and 1.5 mm thick and with 10 or 15 teeth.

3. TECHNICAL SPECIFICATIONS

Code	53001001
Nº gels	1 or 2
Gel dimensions	82x82 mm
Glass plate dimensions	102x85 mm
Nº samples	10 or 15 per gel
Comb thickness	1 and 1.5 mm
Buffer volume (approx.)	500 mL
Working conditions	0-40 °C and < 80% RH
Electrophoresis voltage	200 V
Running time	45-60 min
Cell dimensions (LxWxH)	100x140x150 mm

4. OPERATING MODE

4.1 Preliminary inspection

Unwrap the vertical electrophoresis cell, take off all protections and make sure it does not present any damage because of the shipment. In case it presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- 1 External tank
- 1 Lid with fixed leads
- 1 Gel support
- 2 plastic wedges
- 1 Gel casting device
- 2 Silicon rubber cushions
- Set of glass plates:
 - 4 flat glass plates
 - 2 flat and notched glass plates with 1.0 mm thick fixed spacers
 - 2 flat and notched glass plates with 1.5 mm thick fixed spacers


Components

Nahita vertical electrophoresis cell is supplied together with all the necessary components for preparation and electrophoresis of polyacrylamide gels (82x82 mm):


- *External tank*: made of a unique seamless high quality acrylic material piece to avoid leakage and breakage. It is a durable and transparent device resistant to chemical agents and high pressure and temperatures.

- *Lid*: it is also made of a high quality and resistant material. It presents two fixed power leads with banana connector to be plugged to a power supply.


- *Gel support*: gel casting and electrophoresis are performed at the same support minimizing gel manipulation and speeding up and making work easier. The support includes 2 plastic wedges to hold the gel assembly and pure platinum electrodes resistant to corrosion and high temperatures.

- *Gel casting device*: it is an exclusive device that presents a silicon rubber cushion on its base, allowing a perfect sealing to avoid polyacrylamide leaks that constitute a toxic risk for the user. It also presents two white handles with two positions for casting gels of 1.0 and 1.5 mm thick respectively.

- *Glass plates*: the equipment includes 2 flat glass plates and 2 notched plates with fixed spacers of 1.0 and 1.5 mm thick for casting and electrophoresis of gels. The equipment also includes a thick glass plate to use in case of carrying out a single gel electrophoresis.

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INDEX OF CONTENTS


1. USES OF THE INSTRUMENT	17
2. DESCRIPTION AND COMPONENTS	17
3. TECHNICAL SPECIFICATIONS	19
4. OPERATING MODE	19
5. FORMULATION FOR GELS AND ELECTROPHORESIS BUFFER.....	24
6. TROUBLESHOOTING	29
7. MAINTENANCE AND CLEANING.....	27
ANNEX I: CE CERTIFICATE.....	28

1. USES OF THE INSTRUMENT

Nahita horizontal electrophoresis cell is specially designed for separating proteins and nucleic acids in polyacrylamide gels. The equipment is provided complete with all the necessary accessories for preparation and electrophoresis of 82x82 mm gels.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita vertical electrophoresis cell. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita vertical electrophoresis cell, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the electrophoresis cell avoiding sudden movements, knocks, free fall of heavy / sharp objects on it.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Accessories (including their loss) are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

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